



ICE-PIC 0.5 NTU Calibration/Verification Module

Продукт №:	5222500
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Can be used with the 1720C and D instruments. Verify calibration in less than 60 seconds. No loss of sample flow to complete a verification and no sample cell to recondition after verification.

Saves time

Verification is the fastest available, taking less than one minute.

Ease of use

Compact and lightweight, ready for use at any time.

Accuracy

$\pm 10\%$ or better, and each device is supplied with a Certificate of Accuracy.

Repeatability

Exceptional repeatability performance to be tracked closely.

Affordability

One-time investment, no consumables. One-year warranty.

Спецификации

Какво има в кутията:	ICE-PIC, Instructions.
----------------------	------------------------

Какво има в кутията

ICE-PIC, Instructions.