



Be Right™



ICE-PIC 20 NTU Calibration/Verification Module

Продукт №:	5225000
BGN Цена:	Свържете се с нас
Повикване за дата на доставка	

Verify calibration in less than 60 seconds. Available in 1 and 20 NTU modules. Can be used with the 1720C, D and E instruments. No loss of sample flow to complete a verification. No sample cell to recondition after verification.

Saves time

Verification is the fastest available, taking less than one minute.

Ease of use

Compact and lightweight, ready for use at any time.

Accuracy

±10% or better, and each device is supplied with a Certificate of Accuracy.

Repeatability

Exceptional repeatability performance to be tracked closely.

Affordability

One-time investment, no consumables. One-year warranty.