



## Cal-Cube Calibration/Verification Module

Продукт №:

8220G1300

**ОСТАРЯЛ АРТИКУЛ**

Този артикул вече не е наличен.

It is used to verify the calibration of the ACCU4 and Steady Stream 4 Turbidimeters. Patented glass calibration cube is factory certified to a known USEPA approved formazin standard.

Module verifies the calibration of the Accu4 and Steady Stream 4 Turbidimeters. Patented glass calibration cube is factory certified to a known USEPA approved formazin standard.

---

### Спецификации

Измерване, принцип:

Temporarily insert optional GLI

Cal-Cube™ assembly into sensor and

enter its factory-certified standard value for calibration verification.