



## ComLabo Easy PC SW, за sensION+ GLP, CD, кабел, USB адаптер

Продукт №:

LZW8997.99

**ОСТАРЯЛ АРТИКУЛ**

Този артикул вече не е наличен.

LabCom Easy software, to obtain data of Sension+ GLP instruments from a PC. Kit including CD and cable with 9-contact connector and USB adapter.

- On line graph view of measured value
- Results table personalisation, the user selects the displayed information
- Storage of calibration and measuring results in a data base where the data can be viewed or printed at any moment
- Print report of individual sample or results summary for several samples

Data export to Excel

### **Convenient Data Management.**

Labcom easy allows the acquisition and the management of the obtained results from Sension+ Laboratory instruments.