



TITRAMASTER 85 bidirectional

Продукт №:

R31T022

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Handling sample throughput efficiently while keeping data securely stored for future access is an essential requirement for all laboratories. With TitraMaster 85 user-friendly software, you can program applications directly from the computer while processing data from as many as 7 titration workstations or high-performance ion analysers.

Operator-friendly TitraMaster 85 PC Software is designed to interface the new generation of titration workstations from Radiometer Analytical directly with a computer. It automatically manages application methods and records all data obtained from up to 7