



SAC850 Sample changer, Standard version with cables

Продукт №:

R41T043

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Basic sample changer for automation of small sample series, single-radius sample trays for up to 24 samples.

Single-radius sample trays for up to 24 samples

Parking position on sample tray

Dip rinsing/electrode conditioning performed on 1 or 2 dedicated positions on sample tray

Спецификации

Вид:	Пробосменител
Дължина:	686 mm
Размери (В x Ш x Д):	711 mm x 580 mm x 686 mm