



Be Right™



Hardness (Total & Calcium) Test Kit, Model HA-4P

Продукт №:

145700

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Hardness, Total & Calcium drop count titration/EDTA test kit, range 17.1 - 342 mg/L CaCO₃, 100 tests.

The kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

Utilises drop count titration - a simple, economical titration method

Спецификации

Smallest increment:	17.1 mg/L
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	17,1 - 342 mg/L CaCO ₃
Какво има в кутията:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.
Метод, име:	Drop Count Titration/EDTA
Модел на тестов комплект:	HA-4P
Параметър:	Hardness, Total & Calcium - As CaCO ₃
Платформа:	Drop Count Titration

Какво има в кутията

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.