



Iron Reagent Powder Pills, 0.02-3.00mg/L Fe, FerroVer

Продукт №:
BGN Цена:
Доставка до 2 седмици

85428
Свържете се с нас

Опасен



Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe_3^+) and ferrous (Fe_2^+) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Спецификации

ЕРА compliant:	Да
Бр. тестове:	1000
Измерване, диапазон:	0,02 - 3,00 mg/L Fe
Метод, име:	FerroVer
Метод, номер:	8008
Обем / опаковка, размер:	1000 /pk
Платформа:	Прахови вложки