



DPD Free chlorine, Swiftest dispenser reagent (refill)

Продукт №:

2105560

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



For free chlorine determination by the DPD Free Chlorine method. For use with Hach Swiftest DPD dispenser.

USEPA accepted for reporting for drinking water analysis using Hach Method 8021. Total chlorine and potassium iodide refill vials also available.

Range: 0.02 - 2.00 mg/L Cl₂. Sample size: 10 mL. Approximately 250 tests per vial.

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Спецификации

EPA compliant:	Да
Бр. тестове:	250
Измерване, диапазон:	0,02 - 2,00 mg/L Cl ₂
Метод, име:	DPD
Метод, номер:	8021
Параметър:	Хлор, свободен
Платформа:	Swiftest

Задължителни аксесоари

- Swiftest DPD свободен хлор, диспенсер за реактив и флакон с реактив (Item 2802300)