



EDTA Digital Titrator Cartridge, 0.714 M

Продукт №:

1495901

Опасен

BGN Цена:

Свържете се с нас

Налично



Ready-to-use titrant solution for the analysis of Total Hardness, 10-100+ G.d.H. CaCO₃.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Спецификации

Бр. тестове:	40 - 100
Измерване, диапазон:	10 - 200 G,d,h,
Метод, име:	EDTA (CalVer2)
Метод, номер:	8204
Параметър:	Hardness, Total
Платформа:	Цифров титратор
Разтвор за титруване:	0.714 M EDTA