




Be Right™



Stablcal Turbidity Standards Calibration Kit for 2100P

Продукт №:

2659405

Опасен

BGN Цена:

Свържете се с нас

Налично



Four ready to use, stabilised Formazin Primary standard solutions. Kit includes sealed vials of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

Спецификации

Инструмент:	2100P
Концентрация:	< 0,1, 20, 100, 800 NTU
Обем / опаковка, размер:	Set of 4 sealed vials
Параметър:	Turbidity