



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Продукт №:

6788302

BGN Цена:

Свържете се с нас

Налично

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Спецификации

Какво има в кутията:

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Комплект продукти:

5500 sc Silica Analyser, 2/4 Channels

Какво има в кутията

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions