



CDC565 Conductivity cell, 4-pole, non-platinised, epoxy body, MAB6 plug (Radiometer Analytical)

Продукт №:

E61M003

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Предложени артикули за замяна:

- E61M010

Отговаря на изискванията на USP 24-NF19.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.