



## рНС3085-8 комбиниран рН електрод, температура

Продукт №:

E16M342

**ОСТАРЯЛ АРТИКУЛ**

Този артикул вече не е наличен.

### **The right probe for your application.**

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

### **Glass sensors made for long-life and durable use.**

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.

---

## Спецификации

Гаранция:	6 месеца
Електрод, тип:	рНС3085-8
Материал:	Body: ероху