



Be Right™



Cell adapter, 12 and 13 mm, for 2100N and 2100AN turbidimeters

Продукт №: 3033400

ОСТАРЯЛ АРТИКУЛ
Този артикул вече не е наличен.

Allows use of 12 or 13 mm vials with 2100N / AN laboratory turbidimeters, incl. adapter and vial cover

Used with Hach spectrophotometers and portable laboratories. Makes it possible to use other cells in instruments which normally take a 25-mL cell. For use with Models 2100N and AN.

Спецификации

Кювета, съвместимост:	12 mm and 13 mm vials
Обем на доставката:	Vial adapter, vial cover, instructions
Тегло:	0,05 kg