



Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Продукт №:

224801

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Parameter:	Phosphorus
Smallest increment:	0.02
Бр. тестове:	100
Вид на корпуса:	F
Измерване, диапазон:	0 - 4 mg/L PO ₄
Метод, име:	Color disc/Ascorbic acid
Модел:	PO-19A
Модел на тестов комплект:	PO-19A
	(Includes filtration for turbid samples.)
Параметър:	Phosphorus, Orthophosphate (reactive) - As PO ₄ ³⁻
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄