



Test kit, phosphate, total ortho-/meta-, model PO-24

Продукт №:

225001

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Спецификации

Parameter:	Phosphate
Smallest increment:	0.02
	0.1
	1
Бр. тестове:	100 Clear
	100 Turbid
Вид на корпуса:	F
Измерване, диапазон:	0 - 4 mg/L PO ₄
Метод, име:	Color disc/Ascorbic acid
Модел:	PO-24
Модел на тестов комплект:	PO-24
Параметър:	Phosphate, total ortho-/meta- - As PO ₄ ³⁻
Платформа:	Цвятен диск
Стандартен обхват на измерване:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄