



Iron color disc test kit, model IR-18A

Продукт №:

146500

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Parameter:	Iron
Smallest increment:	0.02
Бр. тестове:	100
Вид:	Две метрични стойности
Вид на корпуса:	B
Измерване, диапазон:	Low Single Range
Метод, име:	Color disc/ 1,10 Phenanthroline
Модел:	IR-18A
Модел на тестов комплект:	IR-18A
Параметър:	Iron, medium range - As Fe**
Платформа:	Цветен диск
Стандартен обхват на измерване:	0.02 - 1.00 mg/L Fe