



Hardness, iron, and pH test kit, model HA-62A

Продукт №:

183701

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Бр. тестове:	100
Измерване, диапазон:	6.6 - 8.4 pH units
Модел:	HA-62A
Параметър:	Hardness, Iron, and pH (6.6-8.4)
Платформа:	Цветен диск
Приложение:	Water Quality