



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Продукт №:

223002

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Спецификации

Бр. тестове:	100
Измерване, диапазон:	0 - 3,4 mg/L
Модел:	CN-39WR
Параметър:	Chlorine, Hardness, Iron, and pH
Платформа:	Цветен диск
Приложение:	Water Quality