



## Test kit, molybdate, model MO-2

Продукт №:

1419301

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

**Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method**

**Effectively monitors corrosion inhibitors**

### Спецификации

Smallest increment:	0.2
	1
Бр. тестове:	100
Вид на корпуса:	F
Измерване, диапазон:	0 - 8 mg/L Mo
Метод, име:	Color disc/ Mercaptoacetic acid
Модел:	MO-2
Модел на тестов комплект:	MO-2
Параметър:	Molybdate - As Molybdenum
Платформа:	Цветен диск