



Be Right™



Test kit, manganese color disc, model MN-PAN

Продукт №:

2350800

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Parameter:	Manganese
Smallest increment:	0.05
Бр. тестове:	50
Вид на корпуса:	D
Измерване, диапазон:	0 - 0,70 mg/L Mn
Метод, име:	Color disc/PAN
Модел:	MN-PAN
Модел на тестов комплект:	MN-PAN
Параметър:	Manganese - As Mn
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 0.70 mg/L Mn