



## Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Продукт №:

143601

Опасен

BGN Цена:

Свържете се с нас

Налично



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO<sub>2</sub>, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

**Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.**

### Спецификации

Smallest increment:	1.25
	2
	5
Бележка под линия:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Бр. тестове:	200
Вид на корпуса:	D
Измерване, диапазон:	5 - 100 mg/L CO <sub>2</sub>
Метод, име:	Drop count titration/  Sodium hydroxide  Phenolphthalein
Модел на тестов комплект:	CA-23
Параметър:	Carbon Dioxide - As CO <sub>2</sub>
Платформа:	Drop Count Titration