



Chloride Low Range Test Kit, Model 8-P

Продукт №:

144001

Опасен

BGN Цена:

Свържете се с нас

Налично



Chloride, single parameter test kit, drop count titration (range: 5 - 100

mg/L) / silver nitrate (range: 20 - 400 mg/L), blue polypropylene case

Using silver nitrate titrant, this kit provides fast, accurate detection of chloride.

Low cost and easy to use drop count titration method using a silver nitrate titrant

Ideal for quickly determining quantitative results

Спецификации

Smallest increment:	5
Бр. тестове:	100
Вид на корпуса:	В
Измерване, диапазон:	5 - 100 mg/L & 20 - 400 mg/L Cl ⁻
Метод, име:	Drop count titration/ Silver nitrate
Модел на тестов комплект:	8-P
Параметър:	Chloride - As Cl ⁻
Платформа:	Drop Count Titration