



## Total Hardness Test Kit, Model 5-EP

Продукт №:

145400

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model 5-EP test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizes drop count titration with a range of 1-30 gpg. Contains approximately 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Спецификации

Smallest increment:	1 gpg
Бр. тестове:	100
Вид на корпуса:	В
Измерване, диапазон:	1 - 30 gpg CaCO <sub>3</sub>
Метод, име:	Drop count titration/ EDTA-powder pillows
Модел на тестов комплект:	5-EP
Параметър:	Hardness, Total - As CaCO <sub>3</sub>
Платформа:	Бюрета