



TEST KIT, 5-EP mg/L HARDNESS

Продукт №:

145401

Опасен

BGN Цена:

Свържете се с нас

Налично



Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	20
Бр. тестове:	100
Вид на корпуса:	В
Измерване, диапазон:	20 - 400 mg/L CaCO ₃
Метод, име:	Drop count titration/ EDTA-powder pillows
Модел на тестов комплект:	5-EP MG-L
Параметър:	Hardness, total - As CaCO ₃
Платформа:	Drop Count Titration