



## Hardness (Total & Calcium) Test Kit, Model HA-4P MG-L

Продукт №:

145701

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Hardness, Total & Calcium (Magnesium by difference) drop count titration/EDTA test kit, range 20 - 400 mg/L CaCO<sub>3</sub>, 100 tests.

The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

**Utilises drop count titration - a simple, economical titration method**

---

### Спецификации

|                           |                                                                                       |
|---------------------------|---------------------------------------------------------------------------------------|
| Smallest increment:       | 20                                                                                    |
| Бр. тестове:              | 100                                                                                   |
| Вид на корпуса:           | D                                                                                     |
| Измерване, диапазон:      | 20 - 400 mg/L CaCO <sub>3</sub>                                                       |
| Какво има в кутията:      | Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case. |
| Метод, име:               | Drop count titration/EDTA                                                             |
| Модел на тестов комплект: | HA-4P                                                                                 |
|                           | MG-L                                                                                  |
| Параметър:                | Hardness, total & Calcium - (Magnesium by difference)                                 |
| Платформа:                | Drop Count Titration                                                                  |

---

### Какво има в кутията

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.