



Free & Total Chlorine Test Kit (Model CN-70)

Продукт №:

1454200

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Parameter:	Chlorine (Free and Total)
Smallest increment:	0.02
	0.1
Бележка под линия:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Бр. тестове:	100 Free, 100 Total
Вид на корпуса:	D
Измерване, диапазон:	0 - 3,4 mg/L Cl_2
Метод, име:	Color disc/DPD
Модел на тестов комплект:	CN-70 ²
Параметър:	Chlorine, low range - As free & total Cl_2
Платформа:	Цветен диск