



Test Kit, chlorine (free and total), model CN-80

Продукт №:

2129000

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Parameter:	Chlorine (Free and Total)
Smallest increment:	0.02
	0.1
	0.5
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	0 - 3,4 mg/L Cl ₂
Метод, име:	Color Disc/DPD
Модел:	CN-80
Модел на тестов комплект:	CN-80 ²
Параметър:	Chlorine, low range - As free & total Cl ₂
	total Cl ₂
	(0-10 mg/L: T only)
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 0.68 mg/L Cl ₂ ; 0 - 3.4 mg/L Cl ₂ ; 0 - 10 mg/L Cl ₂