



Be Right™



Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Продукт №:

2444400

Опасен

BGN Цена:

Свържете се с нас

Налично



Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Спецификации

Smallest increment:	10
Бележка под линия:	*mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Бр. тестове:	100
Вид на корпуса:	B
Измерване, диапазон:	10 - 200 mg/L
Метод, име:	Drop count titration/
	Thiosulfate
Модел на тестов комплект:	CN-21P
Параметър:	Chlorine, mid and high range - As total Cl ₂
Платформа:	Drop Count Titration