



**Be Right™**



## **Dissolved oxygen test kit, model OX-2P**

Продукт №:

BGN Цена:

Не е показана дата на доставка

146900

Свържете се с нас

Опасен



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

## Method / Chemistry: Drop count titration

**Replacement reagent set PN 2439802**

Utilizes drop count titration - a simple, economical titration method

## Спецификации

Smallest increment: 0.2

Бд. тестове:

Вид на корпуса:

Измерване, диапазон: 0.2 - 4.0 mg/L & 1 - 20 mg/L О

Метод: имеем Drop count titration/

### Modified Winkler

## Модел на тестов комплект

0X-2P

## Параметър:

### Dissolved Oxygen - As DO

Платформа: