



AccuVac[®] kit, dissolved oxygen

Продукт №:

2515050

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. The simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Smallest increment:	0.2
Бр. тестове:	25
Вид на корпуса:	G
Измерване, диапазон:	0 - 15 mg/L O ₂
Метод, име:	Color disc/HRDO
Модел:	Color disc/HRDO
Модел на тестов комплект:	AccuVac [®] Kit
Параметър:	Dissolved Oxygen - As DO
Платформа:	Accuvac