



Be Right™



Test kit, sulfate, model SF-1

Продукт №:

225100

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Спецификации

Smallest increment:	50
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	50 - 200 mg/L
Метод, име:	Extinction/Turbidimetric
Модел на тестов комплект:	SF-1
Параметър:	Sulfate - As SO ₄ ²⁻