



Test kit, alkalinity, model AL-DT

Продукт №:

2063700

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Бележка под линия:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Бр. тестове:	100
Вид на корпуса:	F
Измерване, диапазон:	10 - 4,000 mg/L CaCO ₃
Метод, име:	Digital Titrator/Sulfuric acid
Модел на тестов комплект:	AL-DT
Най-малки инкременти, стъпки:	0.1-10
Параметър:	Alkalinity - (P) & Total (MO) as CaCO ₃
Платформа:	Цифров титратор