



**Be Right™**



## Test kit, alkalinity, model AL-TA

Продукт №:

2314500

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

Smallest increment:	23 gpg
Бележка под линия:	385 *mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Бр. тестове:	100
Вид на корпуса:	B
Измерване, диапазон:	23 - 495 gpg & 385 - 8500 mg/L
Метод, име:	Drop count titration/Sulfuric acid
Модел на тестов комплект:	AL-TA
Параметър:	Alkalinity - (P) & Total (MO) As $\text{CaCO}_3$
Платформа:	Drop Count Titration