



Be Right™



Test kit, alkalinity, model AL-AP

Продукт №:

2444300

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	0.4 gpg
Бележка под линия:	1 gpg *mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$, ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Бр. тестове:	100
Вид на корпуса:	B
Измерване, диапазон:	0 - 8 gpg & 1 - 20 gpg
Метод, име:	Drop count titration/Sulfuric acid
Модел на тестов комплект:	AL-AP
Параметър:	Alkalinity - (P) & Total (MO) As CaCO_3
Платформа:	Drop Count Titration