



TPH Immunoassay Reagent Set

Продукт №:

2774300

Опасен

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.



For the semiquantitative determination of total petroleum hydrocarbons (TPH) in soil and water. Please note that this reagent set has changed starting in January, 2017. The antibody cuvettes will now come as 12-mm tubes. For these new tubes, the analyst will need to order a new 8 position tube rack with part number 4879910. If the analyst is utilizing the PCII (pocket colorimeter II) for the TPH test, you will also need a new cell adapter with part number 5954610.

Спецификации

Бр. тестове:	Up to 18
Дълбочина:	3.5 /cm
Измерване, диапазон:	20, 50, 100, 200 ppm threshold
Инструмент:	All except DR3000, DR2000, DR900, DR800 Series
Какво има в кутията:	Includes diluent solution, TPH enzyme conjugate, methanol, color developing solution, stop solution, cuvettes, and four calibrator solutions. For up to 18 tests.
Метод, име:	Semiquantitative Colorimeter Immunoassay
Метод, номер:	10050
Параметър:	TPH (Total Petroleum Hydrocarbons)

Какво има в кутията

Includes diluent solution, TPH enzyme conjugate, methanol, color developing solution, stop solution, cuvettes, and four calibrator solutions. For up to 18 tests.