



Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

Продукт №:

224000

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as NO_2^- and contains 100 tests. Range (mg/L): 0-8/0-1600.

Simple, low-cost colour disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Parameter:	Nitrogen, Nitrite
Smallest increment:	2
	40
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	Single Range
Какво има в кутията:	Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.
Метод, име:	Colour disc/
	Ferrous sulphate
Модел на тестов комплект:	NI-6
Параметър:	Nitrogen, Nitrite, high range - as NO_2^-
Платформа:	Цветен диск

Какво има в кутията

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.