



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Red Cap/ Stem, 6/pk

Продукт №:

2928204

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Спецификации

Stem:	Red Cap/Stem
Еденица:	6 pk
Капацитет:	500 mL
Материал:	LDPE
Обем / опаковка, размер:	6 pk
Цвят:	Clear