



Beaker, polypropylene, low form, 400 mL, 6/pk

Продукт №:

108078

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Спецификации

Диаметър, външен:	99.1 mm with spout
Единица:	6 /pk
Капацитет:	400 mL
Материал:	Polypropylene, Low Form
Материал на тялото (измервателния уред):	Polypropylene