



Be Right™



Beaker, polypropylene, low form, 50 mL, 12/pk

Продукт №: 108071
BGN Цена: Свържете се с нас
Доставка до 2 седмици

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Спецификации

Диаметър, външен:	51.1 mm with spout
Единица:	12 /pk
Капацитет:	50 mL
Материал:	Polypropylene, Low Form,
Материал на тялото (измервателния уред):	Polypropylene