



CEL/890 Advanced Portable Laboratory

Продукт №:

2688500

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Consists of:

- portable colorimeter for chemical water analysis (DR890)
- portable pH/conductivity meter with pH electrode and 4-pole cond. electrode (HQ40D.99.121201),
- digital titrator
- reagents basic equipment
- accessories
- illustrated manuals

Complete with carrying case.

Contains durable instruments, premeasured reagents impervious to environmental contamination, and apparatus and labware

Offers on-site results in an easy-to-use, rugged, portable lab

Rugged kit is designed to withstand harsh field conditions

Hard-sided durable cases

Instruments have option of running on batteries, so tests can be run immediately without concern for sample preservation and handling

versatile for lab or field use

Real-time decisions

In the field, rapid testing and decision making is essential. With this kit, there is no waiting for hours or days for a lab result

Reliable results

Whether used in your mobile lab or at the field site, this CEL890 kit yields consistent results time after time

Спецификации

Какво има в кутията:

The Advanced Portable Laboratory includes DR/890 Colorimeter, HQ40d multi-meter, PHC201 pH electrode (1 meter cable), CDC401conductivity probe (1 meter cable), Digital Titrator, illustrated instrument and procedures manuals, reagent sets, and all required apparatus in a rugged carrying case.

Приложение:

Potable Water

Какво има в кутията

The Advanced Portable Laboratory includes DR/890 Colorimeter, HQ40d multi-meter, PHC201 pH electrode (1 meter cable), CDC401conductivity probe (1 meter cable), Digital Titrator, illustrated instrument and procedures manuals, reagent sets, and all required apparatus in a rugged carrying case.