



Portable Pump

Продукт №:

3206

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

The Sigma Portable Pump uses a peristaltic pump to collect liquid samples from open channels and shallow monitoring wells. It's self-priming and able to run dry with a maximum lift of 27 feet. Housed in a rugged NEMA 4X, 6 fiberglass housing, this pump is

Convenient sample collection

Uses a peristaltic pump to collect liquid samples from open channels and shallow monitoring wells.

Easy to use

Self-priming and able to run dry with a lift of up to 27 feet.

Rugged

Housed in a rugged NEMA 4X, 6 fiberglass housing.

Reliable

Able to withstand temporary submersion.

Спецификации

Какво има в кутията:

Includes: pump, pump tube insert, stainless steel connector with hose clamps, and instructions

Какво има в кутията

Includes: pump, pump tube insert, stainless steel connector with hose clamps, and instructions