



## POCKET Colorimeter II Colorimeter Test Kit for Manganese (HR) analysis

Продукт №: 5870015

Опасен

**ОСТАРЯЛ АРТИКУЛ**  
Този артикул вече не е наличен.



### Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Manganese, 0.20-20.0 mg/L Mn, Periodat Oxidation Method. In carrying case, complete with reagents, cuvettes and operating instructions.

**Simple:** all functions are available via only four keys

**Power:** battery operation for a maximum of 2000 tests

**Clear readings:** even under difficult conditions, thanks to backlit display and large readout

**Reliable results with no mains connection**

**Rugged construction yet lightweight**

---

### Спецификации

|                                              |                                               |
|----------------------------------------------|-----------------------------------------------|
| Smallest increment:                          | 0.1                                           |
| Бр. тестове:                                 | 100                                           |
| Вид на корпуса:                              | D                                             |
| Детектор:                                    | Silicon detector                              |
| Дисплей: течнокристален (фоново осветление): | LCD, backlit                                  |
| Дължина на вълната:                          | Fixed wavelength $\pm$ 2 nm varies with model |
| Измерване, диапазон:                         | 0,2 - 20,0 mg/L Mn                            |
| Метод, име:                                  | Cold periodate                                |
| Метод, номер:                                | 8034                                          |
| Ограждение, клас:                            | IP67 (без отделението за батерии)             |
| Параметър:                                   | Manganese                                     |
| Размери (В x Ш x Д):                         | 155 mm x 61 mm x 35 mm                        |
| Сертификати за съответствие:                 | CE                                            |
| Тегло:                                       | 0,23 kg                                       |