



FerroZine Iron Reagent Solution Pillows, 0.009-1.400mg/L Fe

Продукт №:

230166

Опасен

BGN Цена:

Свържете се с нас

Налично



For the determination of soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron at trace levels.

FerroZine Iron Reagent determines soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron. Digestion is needed for strongly chelated Iron complexes (such as those with EDTA) and Total Iron. For trace level Iron analysis. 50/pack

Спецификации

Бр. тестове:	50
Измерване, диапазон:	0,009 - 1,400 mg/L Fe
Инструмент:	All except DR 900, PC II
Метод, име:	FerroZine
Метод, номер:	8147
Обем / опаковка, размер:	50 /pk
Параметър:	Iron