



Reagent set, copper, porphyrin, 25 mL

Продукт №:

2242700

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Powder reagents in sealed packages for the determination of Copper, Porphyrin Method. For 25ml sample.

For determination of copper by the Porphyrin Method. Hach Method 8143. Range: 2 - 210 µg/L Sample Size: 25 mL. Reagent Set includes: Copper Masking Reagent Powder Pillows, Porphyrin 1 and 2 Reagent Powder Pillows for 100 tests.

Спецификации

Бр. тестове:	100
Измерване, диапазон:	2 - 210 µg/L
Инструмент:	DR/3000, DR/2000
Метод, име:	Porphyrin
Метод, номер:	8143
Платформа:	Прахови вложки