



Copper Masking Reagent Powder Pillows, 10 mL, pk/100

Продукт №:	2603449
BGN Цена:	Свържете се с нас
Налично	

Powder reagent in sealed package for the determination of Copper, Porphyrin Method. For 10ml sample.

For Copper determination by the Porphyrin method. Hach Method 8143 Range: 2 - 210 µg/L Cu. Sample Size: 10 mL Pack of 100 Powder Pillows.

Спецификации

Бр. тестове:	100
Метод, име:	Porphyrin
Обем / опаковка, размер:	/pk
Платформа:	Прахови вложки