



Potassium 3 Reagent Powder Pillows, pk/100

Продукт №:

1432399

Опасен

BGN Цена:

Свържете се с нас

Налично



Powder reagent in sealed package for the analysis of Potassium, 0.1-7.0mg/l K.

For potassium determination in water, soil and plant tissues by the tetraphenylborate method. Hach Method 8049 Sample size: 25 mL. Pack of 100 powder pillows.

Спецификации

Бр. тестове:	100
Метод, име:	Tetraphenylborate
Обем / опаковка, размер:	100 /pk
Платформа:	Прахови вложки