



Be Right™



Sodium Hydroxide Digital Titrator Cartridge, 1.600 N

Продукт №:

1437901

Опасен

BGN Цена:

Свържете се с нас

Налично



Ready-to-use titrant solution for Acid-Base Titration (1-10 meq/l acid) and determination of Acidity (100-4000mg/l CaCO₃).

For acidity determination by Digital Titration, 100 - 4000 mg/L CaCO₃

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Спецификации

EPA compliant:	Да
Бр. тестове:	40 - 100
Измерване, диапазон:	1 - 4000 mg/L CaCO ₃
Метод, име:	Fenolftaleína
Метод, номер:	8233 (Base)
Параметър:	Acidity
Платформа:	Цифров титратор
Разтвор за титруване:	1.600 N NaOH