



Be Right™



Potassium iodide powder pillows for test kits and lab methods, 100 pcs

Продукт №:

107799

Опасен

BGN Цена:

Свържете се с нас

Налично



Powder reagent in sealed packages for the determination of total chlorine with lab titration and chlorine reagent set (drop count).

Reagent for several iodometric titration methods. Pack of 100 powder pillows.

Спецификации

Бр. тестове:	100
Измерване, диапазон:	Various
Метод, номер:	10025
Обем / опаковка, размер:	100 /pk
Платформа:	Прахови вложки